

## GAMARET intramammary suspension

### 1. NAME AND ADDRESS OF THE MARKETING AUTHORISATION HOLDER AND OF THE MANUFACTURING AUTHORISATION HOLDER RESPONSIBLE FOR BATCH RELEASE

Bioveta, a.s.  
Komenského 212,  
683 23 Ivanovice na Hané, Czech Republic

### 2. NAME OF VETERINARY MEDICINAL PRODUCT

GAMARET intramammary suspension

### 3. STATEMENT OF THE ACTIVE INGREDIENTS AND OTHER SUBSTANCES

*1 applicator containing 10 ml of  
homogeneous white to slightly  
yellowish oily suspension contains:*

**Active substances:** Procaini  
benzylpenicillinum monohydricum 100 mg,  
Neomycini sulfas 102 000 IU,  
Dihydrostreptomycini sulfas 125 mg,  
Novobiocinum natrium 100 mg,  
Prednisolonum 10 mg

### 4. INDICATIONS

Treatment of acute and chronic mastitis in  
cows in the period of lactation caused by  
micro-organisms susceptible to novobiocin,  
penicillin, dihydrostreptomycin and  
neomycin.

### 5. CONTRAINDICATIONS

Hypersensitivity to penicillin antibiotics.  
Do not use cleaning wipes on injured  
teats.\*)

### 6. ADVERSE EFFECTS

Mainly after a repeated administration of  
the preparation, there is a possibility in  
hypersensitive individuals to develop  
allergic reactions (salivation, dyspnoea,  
swelling of the skin on the head and  
perineum).

### 7. TARGET ANIMAL SPECIES

Cattle - Milking cows in their lactation  
period.

### 8. DOSAGE FOR EACH SPECIES, ROUTE AND METHOD OF ADMINISTRATION

#### Dosage:

The content of one applicator (10 ml of  
the preparation) is administered to  
the affected quadrant. Before  
the preparation is administered the lacteal  
glands are milking, and thoroughly cleaned  
and disinfected using enclosed cleaning  
wipes.\*)

#### Route and method of administration:

Intramammary. Udder and teats are  
washed as required with warm water and  
thoroughly dried. The ends of teats are  
disinfected with an appropriate agent and  
after removing the plastic cover of  
the applicator's tip is the applicator inserted  
to the teat duct. The content is  
administered to the udder by pressing  
the piston. After the administration is  
the particular quadrant massaged in order  
to ensure a better distribution of  
the preparation into milk cisterns.  
If required by the situation, it is possible to  
repeat therapy after 24 – 48 hours.

### 9. INSTRUCTIONS FOR PROPER ADMINISTRATION

The preparation should be thoroughly  
shaken prior to use!

### 10. WITHDRAWAL PERIOD

Meat – 7 days

Milk – 72 hours (6 x milking).

### 11. SPECIAL PRECAUTIONS FOR STORAGE

Keep out of the reach of children!

Store at a temperate below 25 °C.

Keep in dry place.

Protect from light.

Do not freeze. Do not use after the shelf life  
stated on the cover.

### 12. SPECIAL WARNINGS

Persons handling with the preparation  
should avoid skin or mucous membrane  
contact with the preparation.

In case of an accidental contact, it is  
necessary to wash the site with running  
water.

Individuals with known hypersensitivity to  
penicillin should avoid contact with  
the veterinary medicinal product.

After use cleaning wipes wash your hands.  
For people sensitive or allergic to isopropyl  
alcohol is recommended to use gloves.\*)

### 13. SPECIAL PRECAUTIONS FOR THE DISPOSAL OF UNUSED PRODUCT OR WASTE, IF NECESSARY

All unused veterinary medicinal product or  
waste from this product must be discarded  
according to the local regulations.

### 14. DATE ON WHICH THE PACKAGE INSERT WAS LAST APPROVED

July 2016

### 15. OTHER INFORMATION

For animal treatment only.

This veterinary medical product is available  
only on prescription.

The supplied disinfecting towels are  
moistened with 65% v/v isopropyl alcohol  
solution (2.4 ml/towel).\*)



\*) If the wipes are part of the original packaging

Dimensions (mm):

Colours:

Prepared by:

Date:

A6 (105 × 148 mm) duplex printing

black

Hanáková

24. 6. 2020

Approved by: .....

Date: .....