

LOTAGEN 360 mg/g concentrate for vaginal/cutaneous solution

1. NAME AND ADDRESS OF THE MARKETING AUTHORISATION HOLDER AND OF THE MANUFACTURING AUTHORISATION HOLDER RESPONSIBLE FOR BATCH RELEASE, IF DIFFERENT
Bioveta, a.s., Komenského 212/12,
683 23 Ivanovice na Hané

2. NAME OF THE VETERINARY MEDICINAL PRODUCT
LOTAGEN 360 mg/g concentrate for vaginal/cutaneous solution
Policresulen

3. STATEMENT OF THE ACTIVE SUBSTANCE(S) AND OTHER INGREDIENT(S)
Lotagen is reddish-brown, clear concentrate for vaginal / cutaneous solution, which in 1 g contains:

Active substance:
Policresulen 360 mg/g

4. INDICATION(S)
Sterility

Sterility induced by subacute infections, cervicitis, vaginitis, vulvitis, trichomoniasis.

Obstetrics

Vagina damage, post partum vaginal bleeding, prevention of MMA syndrome (mastitis-metritis-agalactia) in sows.

Surgical interventions and wound treatment

Slight local bleeding spots and bleeding during surgical interventions, fresh and slowly healing wounds, lupus (dermatitis of fetlock in horses), hoof cancer (pododermatitis chronica verrucosa madigans), ulcers on extremities, putrefaction of soft tissues, interdigital necrobacillosis, foot-and-mouth disease in sheep, burns, inflammation of the alvearium (otitis externa), furunculosis, dermal eczema, erosions, irrigation of urovagina in mares, obliteration of lacteal pseudo-fistula.

5. CONTRAINDICATIONS

Not known.

6. ADVERSE REACTIONS

Not known.

7. TARGET SPECIES

Cattle, horses, pigs, sheep, goats, dogs and cats.

8. DOSAGE FOR EACH SPECIES, ROUTE(S) AND METHOD OF ADMINISTRATION

Cattle:

Cervicitis, vaginitis and vulvitis – Irrigate the surface vaginal damages with 2-% solution of the preparation. This procedure is not recommended for deep and perforating wounds.

Lacteal "pseudo-fistula" ("additional lacteal gland") – 9-% solution of the preparation should be applied into a fistula mouth using a milk catheter.

Mares:

Vaginitis (pneumovagina, urovagina) – The vagina should be irrigated with 1-2% solution of the preparation in order to disturb parenchyma. Irrigation of urovagina in horses – 1 - 3 litres of 0.5-% solution of the preparation should be used for irrigation. The procedure can be interrupted only if clear irrigation liquid flows back.

Local Use during Surgical Intervention and Cleaning of Old Wounds:

Hemostasis – In case of slight surface bleeding, cover the wound with gauze soaked with the preparation.

Cleaning of wounds, ulcers, abscesses, eczema and other pathological dermal changes – Apply gauze soaked with 4 – 20-% solution of the preparation onto the old wound or the pathologically changed tissue. Ulcers onto extremities, putrefaction of soft tissues, urovagina, interdigital necrobacillosis, foot-and-mouth disease in sheep, lupus-dermatitis, horn wounds – In case of small surgical interventions, apply 20-% aqueous solution of the preparation.

Use in Small Animals

Eczema on lips and skin folds – Apply gauze soaked with the preparation onto an affected spot.

Repeat the procedure, if necessary.
Fistula of anal gland – Inject 2 ml of 5-% aqueous solution of the preparation into an anal gland cavity.
Repeat the procedure, if necessary.

Interdigital ulceration – Apply gauze soaked with the preparation on the affected spot between fingers.

Otitis externa inflammation – Rinse it with 5-% aqueous solution of the preparation once a day.

The same concentrations (%) and volumes (ml) of the preparation Lotagen conc. sol should be used when the procedures are repeated according to the indications as above.

Preparation of working solutions of Lotagen

Solution concentration (%)	Lotagen (ml)	Purified water
0.5 %	0.5 ml	Up to 100 ml
	5 ml	Up to 1 l
1 %	1 ml	Up to 100 ml
	10 ml	Up to 1 l
2 %	2 ml	Up to 100 ml
	20 ml	Up to 1 l
4 %	4 ml	Up to 100 ml
	40 ml	Up to 1 l
8 %	8 ml	Up to 100 ml
20 %	20 ml	Up to 100 ml

9. ADVICE ON CORRECT ADMINISTRATION

Shake before use.

10. WITHDRAWAL PERIOD

None.

11. SPECIAL STORAGE PRECAUTIONS

Keep out of the reach and sight of children.

Store at the temperature to 25 °C.

Do not refrigerate or freeze.

Do not use after the expiry date stated on the etiquette.

Shelf-life after first opening the container: 21 days.

12. SPECIAL WARNING(S)

After contact with eyes, skin or mucosa may cause irritation. Avoid contact with eyes, skin and mucosa. In case of skin or mucosa affection, immediately wash the site with enough amount of water, remove contaminated clothes and shoes. In case of contact with eyes remove contact lenses and immediately rinse with plenty of water.

Ensure adequate rinsing by manual opening eyes. In case of irritation, seek medical advice. In case of accidental ingestion, rinse mouth with great amount of water, seek medical advice and show the package insert or label to the general practitioner.

13. SPECIAL PRECAUTIONS FOR THE DISPOSAL OF UNUSED PRODUCT OR WASTE MATERIALS, IF ANY

Any unused veterinary medicinal product or waste materials derived from such veterinary medicinal products should be disposed of in accordance with local requirements.

14. DATE ON WHICH THE PACKAGE LEAFLET WAS LAST APPROVED
September 2020



Dimensions (mm): A6 (105 × 148 mm) duplex printing
Colours: black
Prepared by: Hanáková
Date: 01. 10. 2020

Approved by:

Date: